

# Siena

## For Stairs That Turn

Compact and easy to use, the Siena chair from Stannah travels along a rail that is custom-made for your stairs. With over forty years of experience making stairlifts like this, we know how to make the best use of the space available - helping to blend your stairlift discreetly into your life.

Whether you're accessing a basement or making full use of the rooms upstairs, the Siena gives you the freedom to set your own schedule. You can use it to carry the groceries and laundry - and park it neatly out of the way when not in use. Choose from a range of colors, and, if needed, sensible options are available to make this the perfect stairlift for your needs.



**Custom-made** to fit almost any staircase, the Stannah 260 saves you the trouble and cost of moving to an unfamiliar home.



**This versatile design** allows your stairlifts to be placed on either side of the stairs, making the best use of a narrow staircase.



**Unfold** the Siena without bending over, for added safety at both landings.



**Swivel** to face the top landing, or (if space allows) your stairlift can continue on to the top landing.

## Siena 260 Key Features

### Upholstery

Comfortable, waterproof upholstery, available in six colors.

### Chair width options

Standard model is 19½" between the armrests.  
Slim model is 16" between the armrests.

### Easy-to-use controls

The slimline control allows you to use any part of your hand to operate the stairlift.

### Swivel levers

Use the levers, or optional one or two-way powered swivel, to get on and off your stairlift in the safest position at the top or bottom of your stairs.

### Key Lock

The key lock prevents unauthorized use.

### Custom-made

This chair travels along a rail which we custom-make to hug your stairs.



### On/Off switch

Conveniently placed On/Off switch prevents the chair from moving unexpectedly.

### Seat belt

Mounted on the armrest within easy reach. Unfastens with a simple latch.

### Easy folding

The armrests, seat and footrest fold neatly when not in use. The folding mechanism eliminates the need to bend over.

### Safety sensors

These stop the stairlift if it meets an obstruction.

### Armrest safety

The down-side armrest must be lowered when riding the chair.

### Rated load

One person, 300lbs

Model shown: Siena chair with Tan upholstery

## Upholstery Colors



Tan



Ivory



Cranberry



Pine



Walnut



Midnight

# Stannah

### Stannah Stairlifts Inc.

Head Office & Showroom: Tel: (508) 520-7878  
20 Liberty Way, Suite A, Franklin, MA 02038

Tri State Office & Showroom: Tel: (914) 345-8484  
45 Knollwood Rd., Suite 103, Elmsford, NY 10523

E-Mail: [Stannah@stannah.com](mailto:Stannah@stannah.com)

Website: [www.Stannah-Stairlifts.com](http://www.Stannah-Stairlifts.com)